

UDC 57.088 + 57.053

Personalized responsiveness of human PBMCs to the action of IL-7

M. V. Kovalchuk¹, T. A. Ruban¹, M. O. Usenko¹, V. A. Kordium¹¹ Institute of Molecular Biology and Genetics, NAS of Ukraine
150, Akademika Zabolotnoho Str., Kyiv, Ukraine, 03143
kovmv@ukr.net

Interleukin-7 (IL-7) is an important element in the functioning of the immune system. Therefore, it would be appropriate to assess individual responses of the immune cells to IL-7. **Aim.** To analyze *ex vivo* responses of the healthy adult individual T cells to increasing rhIL-7 concentrations. **Methods.** Isolation and cultivation of PBMCs and *in vitro* cytokine bioassay. Quantitative determination of cell viability was performed by the metabolic MTT-uptake assay in viable cells. EC₅₀ was calculated by means of OriginPro7.5. **Results.** The dose-response effect on T cell subsets of PBMCs for different concentrations of rhIL-7 in individuals *in vitro* was evaluated. IL-7 bioassay revealed that the individual EC₅₀ values ranged from 0.41 to 1.5 ng/ml. Inter-individual variability was high for EC₅₀ values (CV_G=38 %). Maximal percentage increase in cell viability among subjects demonstrated moderate variability (CV_G=15 %). **Conclusions.** The study revealed that responsiveness to IL-7 is various for the PBMCs originated from different individuals. Our research assumes that the personalized responsiveness of human PBMCs to IL-7 may be considered as valuable prognostic information about the immune system state and the T cells response to IL-7-based therapies. The advantage of the optimized bioassay for the detection sensitivity of PBMC response to IL-7 was shown.

Key words: IL-7, PBMCs, inter-individual variability.

Introduction

Interleukin-7 (IL-7) is an essential element in the immune system functioning, specifically, in lymphopoiesis and T cell development. The cytokine was initially identified by Namen *et al.* for its ability to stimulate B lymphopoiesis [1]. Now it has been shown that IL-7 plays a critical role in modulating T-cell function *in vivo* [2, 3]. Furthermore, numerous recent re-

searches have demonstrated that IL-7 might serve as an effective adjuvant, which can modulate immune responses and strongly elicit protective immunity [4, 5].

Regulation of the IL-7 action depends on IL-7 receptors (IL-7R) on cells. The IL-7R functioning has some peculiarities. First, whereas the common gamma chain (γ_c) cyto-

kine receptors are up-regulated following T cell activation, the expression of IL-7R alpha chain is lost on the effector cells, but then re-expressed on the memory cells [6, 7]. Second, it has been shown that differences in biological responses may be associated not only with the relative distribution of immune cells, but also with the absolute amount of IL-7R molecules on single cells, which can vary from 18452 to 61228 (mean value) in PBMCs from healthy donors [8]. Third, the amount of IL-7 receptors expressed on a cell is fine-tuned by receptor internalization and suppressor of cytokine signaling proteins. The half-life of IL-7R α decreases from 24 hours to approximately 3 hours after IL-7 treatment [9, 10]. And finally, IL-7 binds strongly to heparin sulfate proteoglycans in the extracellular matrix *in vivo* and does not form a linear gradient from producing cells [11]. Together, these data suggest that IL-7/IL-7R signalling axis participates in many functions of immune system, and perturbation of this pathway is associated with immune disorders.

IL-7 is currently undergoing a number of clinical trials and the use of rhIL-7 has been attempted. The persona characterization of immune cells with respect to cytokine responsiveness and the variation of immune responses among individuals can provide valuable prognostic information for practical use. To date, inter-individual differences in the PBMC responsiveness to IL-7 have not been studied systemically. Thus, the aim of this pilot study was to assess *in vitro* the personalized responsiveness of PBMCs to increasing doses of rhIL-7 and inter-individual variability of the cell responses in healthy individuals.

Materials and Methods

Peripheral blood cell isolation. Cell preparation and culture. Venous peripheral blood samples were obtained from 18 healthy age-matched individuals (age range 24–50 years, males) and processed within 2 h. Primary PBMCs were isolated by Ficoll-Hypaque (Ficoll-Hypaque Plus, GE Healthcare, Sweden) density gradient centrifugation following the manufacturer's instructions and resuspended in complete medium (RPMI 1640 supplemented with 10 % FCS, 50 μ M 2- β -Mercaptoethanol, and penicillin/streptomycin). The viability of mononuclear cells was assessed by the trypan blue dye exclusion test, and it was always more than 95 %.

IL-7 bioassay and assay optimization

The PBMC responsiveness to the action of rhIL-7 was measured in MTT cell viability assay using PHA-activated human peripheral blood lymphocytes. Activated T lymphocytes were obtained by cultivation of PBMCs at a density of 1×10^6 /ml for 3, 5 and 6 days with 10 μ g/ml PHA (Sigma HA-15) in 25-cm² flasks at 37°C 5 % CO₂ in a humidified incubator. Mitogen responses were terminated on days 3, 5, and 6. The PHA-activated cells *in vitro* were harvested and washed three times in cytokine free RPMI + 2 % FCS. 50 μ l of the cell suspension were seeded in triplicates into the wells of the flat-bottomed 96-wells plate (1×10^5 per well). To assess cytokine-mediated effects, the cells were subsequently treated with varying concentrations of rhIL-7 (PeproTech, USA) for 72, 120, and 144 h. The stock treatment IL-7 solution was made up to 1 mg/ml in PBS with 0.1 % BSA. The solution was then adjusted to the indicated concentra-

tion, using culture medium. 50 µl of each treatment solution were then dispensed into each well in triplicate. The final volume in each well was 100 µl and the IL-7 concentration was 62.5, 125, 250, 500, 1000, 2000, 3000, 4000, 5000 and 6000 pg/ml. The plates were taken out after incubation and MTT cell viability assay was carried out [12]. All rhIL-7 samples used in this study were from the same bulk production lot.

Graphs and Statistical Analysis. The dose-response curves were created for each PBMC sample, using the photometric results for different indicated cytokine concentrations. Data were presented as means ± standard deviation (SD). The results of cytokine-induced activation were expressed as cell viability. Quantitative determination of cell viability of PBMCs was performed using metabolic MTT-uptake assay and evaluated comparing cell optical density at different IL-7 concentrations with the same values after starting the cultivation without IL-7 taken as 100 % ($OD_x/OD_{contr} \times 100\%$). Maximal percentage increase in cell viability was expressed as ($OD_{max}/OD_{contr} \times 100\%$). Inter-individual variability was expressed as a coefficient of variation $CV_G(\%) = (1 \text{ SD}/\text{mean}) \times 100\%$ using the overall mean for each variable. The rhIL-7 EC_{50} was defined as the effective concentration of cytokine that elicits a 50 % increase in cell viability in the bioassays. EC_{50} -data were calculated by means of OriginPro 7,5 software using sigmoidal dose-response curves.

Results

A number of assays are focused on exploring effects of the agonists supplement on cellular receptors. The cell functional assays are fre-

quently used to measure cell survival that is a result of downstream effects of the] IL-7R activation. The present study was undertaken to determine whether the individual PBMC responsiveness to rhIL-7, measured as survival of mitogen-stimulated PBMCs *in vitro*, might characterise individual functional status of IL-7 receptors in bloodstream. As a result, some questions emerge. First, what biological variations exist in the PBMC responsiveness among healthy individuals? For that a small-scale age-matched random sample was screened for IL-7 response ability. 18 healthy individuals with normal range of circulating white blood cell counts ($4.0\text{--}10.0$) $\times 10^6$ cells/ml were enrolled in this study.

Second, for noticeable inter-individual difference some modifications of standard IL-7 bioassay were made. Supplementary experiments were designed to establish optimal conditions for getting a standard S-shaped or sigmoid curve used to analyse cell responsiveness to IL-7 since the obtained data on dose-response measurements do not always meet the requirements for the sigmoidality of the curve. We established that the IL-7 dose range 62.5–6000 pg/ml was optimal for the single dose-response curve construction. The IL-7 dose in the range of 0.5–10 ng/ml resulted in the high dose-response curve type but not in S-shaped curve.

Next, the optimal periods between cell collection and PHA response termination as well as between the IL-7 administration and MTT-testing were determined. According to literature, the optimal levels of T cell mitogen responsiveness to PHA in adult PBMCs were after 5 days (optimal on the 6th day) [13]. The time of PHA exposure was selected experi-

mentally. PHA responses were terminated on days 3, 5 and 6 followed by preparation of cells for IL-7 treatment. Using PHA-activated PBMCs which were cultured with PHA for 5 days allowed getting better results for the logistic non-linear curve fit. The time of IL-7 exposure was selected experimentally also. An increase of cell viability was measured in the presence of IL-7 for the same individual PBMCs after 72h and 144-h of IL-7 exposition. The responses of individual PBMCs are shown on representative graphs from 2 experiments (Fig. 1). The EC_{50} values in 144-h culture compared with those in 72-h culture revealed a dose-response shift of 0.943 ± 0.117 to 0.643 ± 0.055 ng/mL for the same PBMCs and the cells became more sensitive to rhIL-7.




Fig. 1. Comparative dose-response effects of rhIL-7 on the identical individual cultured PBMCs *in vitro* following 72-h and 144-h IL-7 exposition period. Cells were stimulated with 0.0625, 0.125, 0.5, 1, 2, 3, 4, 5, and 6 ng/ml of IL-7. Survival data with fitted logistic regression lines are illustrated. Each point represents the mean $(OD_x/OD_{contr}) \times 100\% \pm SD$. A log scale was used for x-axis.

Since the duration of the analysis is important, and the shift in the value of EC_{50} is relatively small, we decided to stay at the 72-h IL-7 exposure. Thus, this assay allowed an assessment of the personalized PBMC responsiveness in the cellular response to IL-7 and evaluation of inter-individual variability.

Personalized cell viability was determined for each individual PBMC specimen and dose-response curves were created. All graphs show the sigmoid curves for rising doses of rhIL-7. The IL-7 concentrations which lead to a 50% increase of cell survival (EC_{50}) were determined for each PBMC specimen and were in the range of 420–1540 pg/ml. The EC_{50} values and maximal percentage increase in cell viability are shown in Table 1.

PBMC responsiveness to IL-7 was heterogeneous between individuals. EC_{50} values demonstrated high inter-individual variability. CV_G was 38%. The difference in sensitivity to IL-7 among individuals was nearly 3 times, based on EC_{50} values. Although a significant variation of PBMC responsiveness to IL-7 was noted, several individual PBMCs were closely approximated to each other according to their functional IL-7 responsiveness that was demonstrated by the overlapping of dose-response curves. Additionally, comparable EC_{50} values were obtained, respectively for specimens N1 and N3 (1.343 ± 0.22 pg/ml and 1.337 ± 0.125 pg/ml), N9 and N17 (767 ± 0.089 pg/ml and 767 ± 0.042 pg/ml), N7 and N15 (1.028 ± 0.213 pg/ml and 1.029 ± 0.21 pg/ml). EC_{50} variability was introduced in Fig. 2 and Table 1.

As for the response magnitude, the percentage increase in cell viability became maximal after 72 h incubation at the IL-7 concentra-

No specimen	EC ₅₀ ±SD (ng/ml)	OD _{max} /OD _{contr} × 100 %
1	1.343±0.22	125
2	1.108±0.123	127
3	1.337±0.125	139
4	1.54±0.2	179
5	0.42±0.09	122
6	1.06±0.10	194
7	1.028±0.213	125
8	0.414±0.17	128
9	0.767±0.089	118

No specimen	EC ₅₀ ±SD (ng/ml)	OD _{max} /OD _{contr} × 100 %
10	0.735±0.096	123
11	0.642±0.071	135
12	0.943±0.117	149
13	0.506±0.106	126
14	0.507±0.123	122
15	1.029±0.21	129
16	1.457±0.163	142
17	0.767±0.042	121
18	1.067±0.073	141

tions 3 ng/ml and more. The IL-7-mediated cell survival of PBMCs was enhanced from 18 % to 94 % in comparison with the untreated cells. Maximal percentage increase in cell viability under IL-7 treatment was not significantly different among individual PBMCs. CV_G for the variable was moderate at 15 %. Heterogeneity in the magnitude of T cell responsiveness to IL-7 among different individuals was presented in Fig. 2 and Table 1 as (OD_{max}/OD_{contr}) × 100 %.

Discussion. IL-7 is currently under extensive investigation as a therapeutic agent in a

variety of immune disorders. In this connection, the biomarkers concerning immune cell responses to IL-7 are needed to monitor the disease and immune status, individual pharmacological responses to IL-7 treatment. The blood is a typical sample obtained in a minimally invasive manner for biomarkers. Nevertheless, cytokines act in either a paracrine or an autocrine manner as they are released and consumed locally, close to the site where the immune reaction occurs. Therefore, they are detectable in peripheral blood at low levels (nearly 2–20 pg in health) [14]. Also,




Fig. 2. Inter-individual variability of the PBMC responses to rhIL-7 *in vitro* (n = 18). (a) Percentages of maximal effects to rhIL-7 in PBMCs (CV_G=15 %). (b) EC₅₀ variability (CV_G=38 %). Each symbol represents the mean±SD of single individual PBMCs.

some studies have demonstrated that there was no significant correlation between serum levels of IL-7 and clinical parameters in patients with several diseases.

As a result, the first question that arises is, “Can the responsiveness of PBMCs to IL-7 be an indicator of the state of the immune system, pharmacological responses to IL-7 treatment on cell physiological level or availability of IL-7 receptors?” It was reported that functional activity of IL-7 receptors can be assessed by the ability of exogenous IL-7 to stimulate the viability (survival + and/or IL-7-induced apoptosis inhibition + proliferation) of freshly isolated immune cells. Such bioassays functionally assess the immediate cellular response to cytokine. To our knowledge, activation of T lymphocytes occurs in lymphoid tissue and generates functionally diverse subsets of effector T cells with the capacity to migrate to multiple tissue sites.

Unfortunately, in human, an assessment of immune cells in tissue sites (lymphoid or injured) is limited. Consequently, the next question that arises is “What types of circulation subsets of immune cells could respond to exogenous IL-7 in bioassay *in vitro*?” In general, an analysis of the blood CD4⁺ T and CD8⁺ T cells subsets as major responders to IL-7 revealed that they contained naïve, central memory (T_{cm}), effector memory (T_{em}), and terminally differentiated effector of T cell populations (T_{emra}) [15].

It should nevertheless be noted that several peripheral blood T cell subsets remain unresponsive to IL-7 stimulation alone throughout the culture period [2, 16]. For these reasons, some ways exist to study cell responsiveness to IL-7 using T cells. Some *in vitro* systems

could reproduce the *in vivo* models illustrating IL-7-related modulation of peripheral T cell homeostasis or even immune processes in local sites with high IL-7 levels. According to the literature two main types of study design prevail to assess T cell responsiveness to IL-7. On the one hand, because homeostatic survival and proliferation of T cells are dependent on signalling from TCR as well as on the IL-7R activation by IL-7, the so called one-step costimulation assays with PHA plus IL-7 are often presented in literature. Human cells were co-stimulated with suboptimal concentrations of mitogen Leucoagglutinin PHA-L to enable proliferative responses to IL-7 [2, 17].

On the other hand, another stimulation scheme exists to assess IL-7 effect. It is a step-by-step approach. The first step is the PHA-driven culture. Only a small part of peripheral blood lymphocytes spontaneously divide *in vitro* but the majority enter the cell cycle following lectin stimulation. The PHA-activated expanded T cells are washed by cytokine-free medium and used for the second step of the assay for IL-7-driven culture. Then the cells are re-cultured in medium contained rhIL-7. We chose the second stimulation scheme, which allows us to use simultaneously the cells of the first step of the assay for other assays. Two parameters were found to be critical for detection of sensitivity in two-step assay and for inter-individual discrimination: prolongation of the time of the PHA-culture; and using low IL-7 doses (62.5–6000 pg/ml). There is an experimental evidence that after 5 days of PHA treatment some metabolic events lead to contraction of expanded T blasts. The expansion of activated T cells is accompanied by an increased sensitivity to apoptosis, which regu-

lates proliferation and maintains lymphocyte homeostasis [13, 18, 19]. On the basis of our results we also hypothesize that some activated T cells after 5 days of PHA-driven culture would be maintained by exogenous IL-7 due to apoptosis regulation. The same cells would be also more sensitive to low concentrations of IL-7. High doses of IL-7 *in vitro* have been demonstrated to stimulate T cell proliferation whereas low doses - to maintain homeostatic T cell survival. Noteworthy, rhIL-7 concentrations above 6 ng/ml did not improve the survival of the tested PBMCs. Moreover, in some cases, the cell survival at 5–6 ng/ml was lower than at 4 ng/ml and lower values of IL-7 concentration. This could be due to the fact that prolonged cultivation of T cells at high IL-7 concentrations results in the internalization of IL-7R on T cells.

There is significant heterogeneity in the immune responses in humans but the extent of variation is only starting to become clear [20, 21]. We tried to estimate inter-individual variability in PBMC responses to IL-7. The responses tended to vary markedly across individuals. There was substantial difference between healthy human subjects in the sensitivity to IL-7 according to values EC_{50} . CV_G for the variance was 38 %. Furthermore, high inter-individual variability is characteristic for other parameters of the immune system. For example, based on the desirable biological variation database specifications www.westgard.com/biodatabase1.htm, the variation coefficients (CV_G) of 35.3 %, 76.4 %, 47.3 % were reported for the lymphocyte count, eosinophils count and immunoglobulins M, respectively. The magnitude of the response, expressed as the percentage increase in cell viability, had

less inter-individual variation ($CV_G = 15 \%$) and was more appropriate for practical use.

Thus, the advantage of the optimized stimulation scheme on the detection sensitivity of PBMC responsiveness to rhIL-7 was shown. The present study revealed that the responsiveness to IL-7 is different for PBMCs originated from various individuals. The research suggests that the personalized responsiveness of human PBMCs to IL-7 may be considered as valuable prognostic information about the immune system state and treatment response in T cell subsets to IL-7-based therapies after confirmation by further validation in large-scale study.

REFERENCES

1. Namen AE, Lupton S, Hjerrild K, Wignall J, Mochizuki DY, Schmierer A, Mosley B, March CJ, Urdal D, Gillis S. Stimulation of B-cell progenitors by cloned murine interleukin-7. *Nature*. 1988; **333**:571–573.
2. Churchman SM, El-Jawhari JJ, Burska AN, Parmar R, Goëb V, Conaghan PG, Emery P, Ponchel F. Modulation of peripheral T-cell function by interleukin-7 in rheumatoid arthritis. *Arthritis Res Ther*. 2014; **16**(6):511–524.
3. Kondrack RM, Harbertson J, Tan JT, McBreen ME, Surh CD, Bradley LM. Interleukin 7 Regulates the Survival and Generation of Memory CD4 Cells. *J Exp Med*. 2003; **198**(12):1797–1806.
4. A Reum Kim, Junsik Park, Jong Hoon Kim, Jeong-Eun Kwak, Youngran Cho, Hyojin Lee, Moon-sup Jeong, Su-Hyung Park, Eui-Cheol Shin. Herpes Zoster DNA Vaccines with IL-7 and IL-33 Molecular Adjuvants Elicit Protective T Cell Immunity. *Immune Netw*. 2018. **18**(5): e38.
5. Nam HJ, Song MY, Choi DH, Yang SH, Jin HT, Sung YC. Marked enhancement of antigen-specific T-cell responses by IL-7-fused nonlytic, but not lytic, Fc as a genetic adjuvant. *Eur. J. Immunol*. 2010; **40**(2):351–358.

6. Jianbao Gao, Lintao Zhao, Yisong Y Wan, Bo Zhu. Mechanism of Action of IL-7 and Its Potential Applications and Limitations in Cancer Immunotherapy. *Int J Mol Sci*. 2015; **16**(5):10267–10280.
7. Mazzucchelli R, Durum SK. Interleukin-7 receptor expression: intelligent design. *Nat Rev Immunol*. 2007; **7**:144–154.
8. Vudattu NK, Kuhlmann-Berenzon S, Khademi M, Seyfert V, Olsson T, Maeurer MJ. Increased Numbers of IL-7 Receptor Molecules on CD4+CD25–CD107a+ T-Cells in Patients with Autoimmune Diseases Affecting the Central Nervous System. *PLoS ONE*. 2009; **4**(8):e6534. doi:10.1371/journal.pone.0006534.
9. Hengires CM, Rino J, Nibbs RJ, Graham GJ, Barata JT. IL-7 induces rapid clathrin-mediated internalization and JAK3-dependent degradation of IL-7R α in T cells. *Blood*. 2010; **115**(16):3269–3277.
10. Faller EM, Ghazawi FM, Cavar M, MacPherson PA. IL-7 induces clathrin mediated endocytosis of CD127 and subsequent degradation by the proteasome in primary human CD8 T cells. *Immunol. Cell Biol*. 2016; **94**:196–207.
11. Kang J, Coles M. IL-7: the global builder of the innate lymphoid network and beyond, one niche at a time. *Semin Immunol*. 2012; **24**(3):190–197.
12. Van de Loosdrecht AA, Beelen RH, Ossenkoppele GJ, Broekhoven MG, Langenhuijsen MM. A tetrazolium-based colorimetric MTT assay to quantitate human monocyte mediated cytotoxicity against leukemic cells from cell lines and patients with acute myeloid leukaemia. *J Immunol. Methods*. 1994; **174**(1–2):311–320.
13. Griffin JF, Bulmer R, Wilson EW. Cord blood lymphocyte subpopulations and mitogenic activity in whole blood microculture. *J Reprod Immunol*. 1979; **1**(4):219–27.
14. Burska A, Boissinot M, Ponchel F. Cytokines as Biomarkers in Rheumatoid Arthritis. *Mediators of Inflammation*. 2014; Article ID 545493, 24 pages <http://dx.doi.org/10.1155/2014/545493>.
15. Sathaliyawala T, Kubota M, Yudanin N, Turner D, Camp P, Thome JC, Bickham KL, Lerner H, Goldstein M, Sykes M, Kato T, Farber DL. Distribution and compartmentalization of human circulating and tissue-resident memory T cell subsets. *Immunity*. 2013; **38**(1):187–197.
16. Hassan J, Reen DJ. Human Recent Thymic Emigrants—Identification, Expansion, And Survival Characteristics. *Journal Immunol*. 2001; **167**(4):1970–1976.
17. Cote S, Matte J, Sad S, Angel JB, Crawley AM. Complexed soluble IL-7 receptor α and IL-7 increase IL-7-mediated proliferation and viability of CD8+ T-cells *in vitro*. *Cell Immunol*. 2015; **293**(2):122–125.
18. Akbar AN, Salmon M. Cellular environments and apoptosis: tissue microenvironments control activated T-cell death. *Immunol. Today*. 1997; **18**:72.
19. Bortwick NJ, Lowdell M, Salmon M, Akbar AN. Loss of CD28 expression on CD8+ T cells is induced by IL-2 receptor γ chain signalling cytokines and type I IFN, and increases susceptibility to activation-induced apoptosis. *Int Immunol*. 2000; **12** (7):1005–1013.
20. Brodin P, Davis MM. Human immune system variation. *Nat Rev Immunol*. 2017 Jan; **17**(1):21–29.
21. Najib Aziz, Roger Detels, Joshua J. Quint, David Gjertson, Timothy Ryner, Anthony W. Butch. Biological variation of immunological blood biomarkers in healthy individuals and quality goals for biomarker tests. *BMC Immunology*. 2019; **20**:33 doi.org/10.1186/s12865-019-0313–0.

Персоналізована відповідь мононуклеарних клітин периферичної крові людини (МНПК) на дію ІЛ-7

М. В. Ковальчук, Т. П. Рубан, М. О. Усенко, В. А. Кордюм

Інтерлейкін-7 (ІЛ-7) є важливим елементом функціонування імунної системи. Тому було б доцільно оцінити індивідуальні реакції імунних клітин на ІЛ-7. **Мета.** Проаналізувати *ex vivo* реакції Т-клітин МНПК здорових дорослих індивідів до зростання концентрації rhIL-7. **Методи.** Виділення та культивування МНПК та аналіз біологічної активності цитокінів *in vitro*. Кількісне визначення життєздатності клітин проводили за допомогою МТТ-тесту у живих клітинах. EC_{50} розраховували використовуючи програмне забезпечен-

ня OriginPro7,5. **Результати.** Було кількісно оцінено вплив різних концентрацій rhIL-7 на Т-клітини МНПК *in vitro* у різних індивідів. Аналіз біологічної активності ІЛ-7 виявив, що індивідуальні значення EC_{50} знаходилися в межах 0,41–1,5 нг/мл. Міжіндивідуальна мінливість була високою для значень EC_{50} ІЛ-7 ($CV_G=38\%$). Максимальний відсоток підвищення життєздатності клітин серед суб'єктів демонстрував помірну мінливість ($CV_G=15\%$). **Висновки.** Дослідження показало, що існує різниця у відповіді МНПК, отриманих від різних індивідів, до дії ІЛ-7. Наші дослідження припускають, що персоналізована відповідь МНПК людини до ІЛ-7 може вважатися корисною прогностичною інформацією про стан імунної системи та реакцію Т-клітин на лікування на основі препаратів ІЛ-7. Показано перевагу оптимізованого методу за чутливістю детектування індивідуальної відповіді МНПК до ІЛ-7.

Ключові слова: ІЛ-7, МНПК, міжіндивідуальна мінливість.

Персонализированный ответ мононуклеарных клеток периферической крови человека (МНПК) на действие ИЛ-7

М. В. Ковальчук, Т. А. Рубан, М. А. Усенко,
В. А. Кордюм

Интерлейкин-7 (ИЛ-7) является важным элементом функционирования иммунной системы. Поэтому было бы целесообразно оценить индивидуальные реакции иммунных клеток на ИЛ-7. **Цель.** Проанализировать

ex vivo реакции Т-клеток МНПК здоровых взрослых индивидов к возрастанию концентрации rhIL-7. **Методы.** Выделение и культивирование МНПК и анализ биологической активности цитокинов *in vitro*. Количественное определение жизнеспособности клеток проводили с помощью МТТ-теста в живых клетках. EC_{50} рассчитывали используя программное обеспечение OriginPro7,5. **Результаты.** Было количественно оценено влияние различных концентраций rhIL-7 на Т-клетки МНПК *in vitro* у разных индивидов. Анализ биологической активности ИЛ-7 обнаружил, что индивидуальные значения EC_{50} находились в пределах 0,41–1,5 нг/мл. Межиндивидуальная изменчивость была высокой для значений EC_{50} ИЛ-7 ($CV_G=38\%$). Максимальный процент повышения жизнеспособности клеток среди субъектов демонстрировал умеренную изменчивость ($CV_G=15\%$). **Выводы.** Исследование показало, что существует разница в ответе МНПК, полученных от разных индивидов, к действию ИЛ-7. Наши исследования предполагают, что персонализированный ответ МНПК человека к действию ИЛ-7 может считаться полезной прогностической информацией о состоянии иммунной системы и реакции Т-клеток на лечение, основанное на препаратах ИЛ-7. Показано преимущество оптимизированного метода по чувствительности детектирования индивидуального ответа МНПК к ИЛ-7.

Ключевые слова: ИЛ-7, МНПК, межиндивидуальная изменчивость.

Received 28.04.2020